

# (Particle sizes for Danish sediments)

Sediment particle size data from 6 size ranges at 3 sites (A, B, C), at 2 depths (2m and 5m), and 5 replicate samples from each site/depth combination are in C:\Examples v7\Denmark PSA\ **Denmark PSA histogram**. The data are already standardised to % composition – if not, you would need **Pre-treatment>Standardise>**(Standardise•Samples)&(By•Total). So, take **Pre-treatment>Cumulate Samples>**(Variable Order•As worksheet). The variable labels are now inaccurate, so replace them by copying/pasting **cumulative** from **Edit>Indicators** and into the **Edit>Labels>Variables** dialog.

---

Revision #2

Created 22 May 2024 01:25:50 by Arden

Updated 22 May 2024 01:29:28 by Arden